



## QUOTATION

To: SAEON Elwandle Node  
ATTN: DR SHAUN HP DEYZEL  
TEL No:  
EMAIL: shaun@saeon.ac.za

DATE: 29 July 2013  
NO OF PAGES: 3  
REF: PE 2995/13 saeon  
FROM: Quanitah Odayar

Dear Shaun,

Thank you for your enquiry. We herewith take great pleasure in quoting as follows:

### OLYMPUS MICROSCOPE

#### SZX10

##### Basic Configuration with Reflected & Transmitted Light

N1977200	SZX2-ZB10	Stereo zoom body with Gallileo type optical system. Zoom ratio 10 : 1, zoom range 6.3 - 63 (1x objective, 10x eyepiece), field number 22. Build in zoom click stop with 11 positions at 0.63/ 0.8/ 1/ 1.25/ 1.6/ 2/ 2.5/ 3.2/ 4/ 5/ 6.3, aperture stop incorporated.
N1978900	SZX2-FOF	Manual focussing unit with coaxial coarse and fine focus for SZX2 zoom bodies. Mount for 32mm Ø pillars (SZX2-ILL#, SZX2-ST#). Focus stroke 80mm, with 36.8mm stroke per rotation for coarse focus and 0.77mm per rotation for fine focus. Maximum load capacity 20kg.
33714	SZX-BI30	Erect image binocular tube with inclination of 30° for SZX7/10 stereo microscopes. Maximum field number 22 (compatible with WHN and WHSZ eyepieces). Interpupillary distance adjustment 51-76 mm (Siedentopf type).
N1198000	WHSZ10X-H	Eyepiece 10x with ESD capability, FN22, focusable, with Diopter adjustment -8 to +5, accepts eyepiece reticule: 24mm 1.5mm thickness
N1200500	DFPLAPO1X-4	Distortion free plan apochromat stereo objective for SZX7 and SZX10. Magnification 1x, numerical aperture 0.1, working distance 81 mm.
N2211800	SZ2-ILST-5	Transmitted and reflected (oblique) light LED illuminator stand for brightfield observation with SZ61/51 and SZX7, adjustable reflected light inclination angle,

#### **Wirsam Scientific & Precision Equipment (Pty) Ltd.**

1 King Edward Street, Newton Park, Port Elizabeth 6045. P.O. Box 27929 Greenacres 6057.

Tel: +27 (0) 41 365 2060 Fax: +27 (0) 41 365 2143

e-mail: [wirsampe@wirsam.com](mailto:wirsampe@wirsam.com) web: [www.wirsam.com](http://www.wirsam.com)

Director: PF Wirsam

	with glass stage plate, color temperature approx. 7000k, focussing stroke 120mm, to adapt SZ61/51 or SZX7 zoom bodies, load capacity of the focusing arm: 7 kg
N2190500 U-LS30-5	30 W Halogen lamphouse for BX43F, BX46F microscope frames and SZX2 ILLK, SZX2-ILLB illumination stands. U-LS30ADP adapter and U-TL4-5 power supply is required for adaptation on BX3 microscope frames.
E0430150 6V30W	6V30 W Halogen Bulb
35425 UYCP	Power Cord
35440 BO61 DUST COV BH2/SZH	Dust cover for BHS/T, BHSU, BHTU, BHSP, BHSM, SZH, BHM, BHA/B/C, CX, SZX12/9

**TOTAL R 54 394.11**

## **Terms and Conditions**

### **VALIDITY**

The prices are valid for 30 days from the date of the quotation and thereafter subject to written confirmation. Any alteration to this quotation shall not be binding unless agreed to in writing by the seller. Prices quoted are subject to fluctuation in the foreign rate of exchange.

The following rates of exchange were applied: Euro1.00 = ZAR13.22. Any variation between the rate quoted and the rate paid to our suppliers will be for your account. Changes in statutory duties and surcharges will also be for your account. In this instance the difference in price will be reflected on your Invoice.

### **Prices**

Our quoted prices **do not** include vat.

### **Delivery**

Approx 4-6 weeks from receipt of official order.

### **Guarantee**

Equipment offered in this quotation is guaranteed against defective material and workmanship from the date of installation for a period of 12 months, except wear and tear parts and parts subject to consumption. Any labour under guarantee will be carried out during normal working hours.

### **Payment**

Prepayment on ordering unless account is held.

### **Cancellation of orders**

Orders are firm and not subject to cancellation. Where cancellation of an order is accepted the purchaser will become liable for a cancellation fee of not less than 20% (twenty percent) of the selling price, plus the cost of shipping the goods back to the factory in the Country of Origin.

### **Installation**

Prices quoted exclude installation.

### **Service**

Complete service facilities are available from our Head Office in Johannesburg as well as our Branch offices. Service contracts are also available.

Wirsam Scientifics' standard Terms and Conditions apply, which are available on request.

---

Should you require any further information, please do not hesitate to contact us.  
We look forward to hearing from you in the near future.

Yours sincerely,

**Quanitah Odayar**

**Wirsam Scientific & Precision Equipment (Pty) Ltd.**

1 King Edward Street, Newton Park, Port Elizabeth 6045. P.O. Box 27929 Greenacres 6057.

Tel: +27 (0) 41 365 2060 Fax: +27 (0) 41 365 2143

e-mail: [wirsampe@wirsam.com](mailto:wirsampe@wirsam.com) web: [www.wirsam.com](http://www.wirsam.com)

Director: PF Wirsam